

DRIGALSKI LACTOSE AGAR

A lactose containing selective and differential medium for detection and enumeration of coliforms.

Dehydrated media	
Code number:	500 g: DAS20500, 5 kg: DAS25000
Colour:	Beige
Appearance:	Homogeneous hygroscopic powder
pH before sterilization (25 °C):	7,1 - 7,5

Direction: Suspend **66 g** in one litre of distilled water and heat with frequent agitation until the medium boils well. Mix well before pouring.

Warning!

The medium is heat sensitive.
No further sterilisation is necessary or desirable.

Prepared media	
Bottled media:	100 ml: DAS30100, 500 ml: DAS30500
Plated media:	55 mm: DAS50055, 90 mm: DAS50090
Colour:	Green
pH (25 °C):	7,2 - 7,4

Direction: Dispense the melted bottled media aseptically into sterile Petri-dishes. Media in Petri-dishes are ready to use.

FORMULA in g/l

Peptones	17,400
Bile salts	1,500
Sucrose	17,000
Lactose	13,000
Sodium chloride	2,000
Bromothymol blue	0,150
Crystal violet	0,002
Agar	15,000

Note: The typical formula can be adjusted to obtain optimal performance.

Storage conditions: Store the dehydrated media tightly closed in a dry place at room temperature. Store the bottled media protected from light at room temperature. Store the plated media protected from light at 2-8 °C. Use before the expiry date on the label.

Quality control:

Test strains	Incubation temp: 37 °C	Growth	Incubation time: 24 h
<i>Escherichia coli</i>	ATCC 25922	Good, yellow colonies	
<i>Salmonella typhimurium</i>	ATCC 14028	Good, blue colonies	
<i>Enterococcus faecalis</i>	ATCC 29212	Inhibited	

References: Ewing (1986) Edwards and Ewing's identifications of the enterobacteriaceae, 4th ed.

In vitro diagnostic – for professional use only!